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## Instruction For Use

valid for all products **Pleura Puncture Set (non-sterile/BULK)**

1. Manufacturer: MEDINORM Medizintechnik GmbH, Gewerbepark 7-9, 66583 Spiesen-Elversberg, Germany
2. **Pleura Puncture Set**; for acute use;
3. Intended use of sterile device: Secretion derive from Cavitas Puralis. Indication: Discharge of Cavitas Pleuralis, Specimen. Contraindication: disordered blood clotting, Pleura Abszess. Specification of patient target group only regarding to Medical Result.
4. **Non-Sterile Medical Device**; Unpack under Cleanroom condition; Do not use by damaged inner packaging! For single use in sterile version!
5. Residual risk or risk by unintended use: risk of infection!
6. Protect from strong light!
7. Unused product is recyclable; used product is to be disposed contamination-protected!
8. Risk of infection! Risk of stitch damage! Observe National Pinprick Guidelines and Regulations!
9. In case of incidence consider national notification requirements!
10. Use through professional persons only

**Note:** Intended use only possible STERILE!

**Note:** Necessity of existing validated sterile packaging process!

**Note:** Sterilisation (EO) through validated process necessary!

**Note:** Sterilisation once possible under use of below instructed parameters!

**Note:** Repeated or Re-Sterilisation must be validated separately!

### Process parameter to be used:

<b>Preconditioning:</b>	
Chamber temperature at end of conditioning phase(°C)	40 ±5
Relative humidity in chamber end of conditioning phase (rH%)	60 ±15
Preconditioning time (h)	6 – 48
<b>Sterilisation:</b>	
Initial vacuum (mbar <sub>abs</sub> )	100 +0/-30
Relative humidity in chamber (at the end of conditioning) (rH%)	70 ±20
Pressure increase after EO injection (mbar <sub>abs</sub> )	390 +30/-0
Gas exposure time (min)	240 +5/-0
Chamber temperature during gas exposure time (°C)	40 +6/-0
<b>Post conditioning:</b>	
Temperature post conditioning chamber (°C)	39 ± 6
Time of desorption (h)	≥48

Explanation of symbols used: <https://medinorm.eu/IFU/Symbols.pdf>

This user manual can be requested in paper form from the manufacturer MEDINORM free of charge. Delivery takes place within seven calendar days. Always have a current version in an appropriate form in your emergency management system!